


NOW IS THE TIME TO REFRESH YOUR OFFICE FOR A NEW START IN 2025

Hausted
by graham-field

TREATMENT TABLES

Hausted treatment tables are designed to provide a comfortable, and supportive base for treatment and diagnostics in an array of medical settings. These tables are available in static or adjustable models with or without pass through drawers.



4100
SERIES
flat top
*4100-xx shown in aqua green | AQ



4103
SERIES
shelf, 2 pass-through drawers, and flat top
*4103-DD-xx shown in dusty blue | DB



4105
SERIES
shelf, adjustable back and contour top
*4105-xx shown in pearl | PL
shelf, 2 pass-through drawers, adjustable back and contour top
*4105-DD-xx†

Hausted
by graham-field

EXAM ROOM TREATMENT TABLE

All 4200 series exam tables have three metal pass-through drawers rather than the typical plastic, so items can be accessed on both sides

- Two Storage and Three Pass Through Drawers
- Pull-Out Leg Rest with Integrated Treatment Tray
- Five Position Continuous Lock Stirrups
- All-metal base table construction
- Powder Coat Finish - Foghorn
- Pelvic Tilt



4200
SERIES
4200-xx†
shown in saddle tan | ST
Options: 4200-E0-xx Integrated Outlet
*† 4200-DW-xx Drawer Warmer
+ Integrated Outlet

Hausted
by graham-field

HIGH / LOW POWER EXAM TABLE

The medical exam table is essential in any health care facility. Hausted provides high-quality options to meet the needs of hospitals, clinics, ERs and more.

- An all-metal construction for the base table provides longevity
- Pull-Out Leg Rest with Integrated Treatment Tray
 - ◊ Tray is independent of the Legrest for any time, easy access
- Hand Pendant controls hi-low, back, flat, return to chair
- UL listed (60601-1) DC electromechanical drive
- Five Position Continuous Lock Stirrups
- Powder Coat Finish - Foghorn
- Power Assisted Back
- Pelvic Tilt

4600
SERIES

*4600-xx
shown in pearl | PL
and in artichoke | AR


Options: 4600-FS-xx Foot Switch
*† 4600-AA-xx Assist Arms
4600-PS-xx Storage Drawer
4600-PS-DW-xx Drawer with Warmer
4600-FL-xx Fully Loaded with All Options§



Hausted
by graham-field

REPLACEMENT TOPS

Make your exam tables look new at a fraction of the cost. Our replacement tops fit many major brand tables.



100
*GF100RT-xx
shown in violet | VI

104
*GF104RT-xx
shown in teal | TE

204
*GF204RT-xx
shown in british tan | BT

304
*GF304RT-xx
shown in ocean gray | OG

404
*GF404RT-xx
shown in doe | DO

405
*GF405RT-xx
shown in shadow | SH

420
*GF420RT-xx
shown in artichoke | AR

460
*GF460RT-xx
shown in british tan | BT

EXAM ROOM FURNISHINGS

A COMPLETE EXAM ROOM SOLUTION

Lumex® products complete the Hausted® Exam Room Furnishings line and are available in color matched vinyl upholstery. DEHP-free vinyl covering available in 24 industry standard colors. See our online swatch card here [📄](#)

LUMEX®
 by graham-field **BARIATRIC PHLEBOTOMY CHAIRS**

The Lumex 6200 Series bariatric blood draw/phlebotomy chair is designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

- Fully upholstered padded seat, back, and armrests
- Powder coat frame finish - Foghorn
- Individual height adjustable arm rests
- Single Flip Arm standard
- Dual Flip Arm option available
- Shown with wall saver wireframe accessory[§] (6200-WALL SAVER)²



6200
SERIES
 Single Flip Arm
 shown in british tan | BT
 ‡6200-xx

800
 POUND
 WEIGHT CAPACITY
 EVENLY DISTRIBUTED



6200
SERIES
 Dual Flip Arm
 shown in blue | BU
 ‡6200DA-xx†

LUMEX®
 by graham-field **PHYSICIAN STOOLS**

We offer a variety of medical stools to meet the needs of an active facility. Offering preferred stool features manufactured in the USA with high quality materials to ensure a long life.

9000 SERIES
 Physician Exam Stool

- 16 inch (40.64 cm) seat diameter
- Choice of three base styles
- Choice of three hand controlled pneumatic assist height adjustment

7000 SERIES
 Physician Surgical Stool

- 16 inch (40.64 cm) seat diameter
- Brushed aluminum base
- Foot-controlled hydraulic height adjustment

350
 POUND
 WEIGHT CAPACITY
 EVENLY DISTRIBUTED



9501
SERIES
 Black Composite Base
 5 configurations available
 shown in gray | GR
 ‡9501-xx



9702
SERIES
 Toe Cap Aluminum Base
 5 configurations available
 shown in plum | PU
 ‡9702-xx



9903
SERIES
 Polished Aluminum Base
 3 configurations available
 shown in british tan | BT
 ‡9903-xx



7408
SERIES
 Brushed Aluminum Base
 3 configurations available
 shown in pearl | PL
 ‡7408-xx

Options: ★†

with Backrest **9501-B-xx**
 with Footring **9501-FR-xx**
 with Backrest & Footring **9501-BFR-xx**
 Lab Height **9501-10BFR-xx**
 with Backrest and Footring

with Backrest **9702-B-xx**
 with Footring **9702-FR-xx**
 with Backrest & Footring **9702-BFR-xx**
 Lab Height **9702-10BFR-xx**
 with Backrest and Footring

with Footring **9903-FR-xx**
 Lab Height with Footring **9903-10FR-xx**

with Backrest **7408-B-xx**
 with Crescent Backrest **7408-CB-xx**



GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

•1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
 Hausted is a registered trademark of GF Health Products, Inc.

Hausted®
 by graham-field **SIDE CHAIRS**

Looking for comfortable and contemporary patient seating for your lobby or waiting room? Our wall saver design keeps the backrest further away from the wall to reduce scratching or marring keeping maintenance costs down and your facility looking its best.

- 14 Gauge Steel Tubing
- Powder Coat Finish-Black
- Wall Saver Leg Design

2000
SERIES
 Wall Saver Side Chair
 shown in violet | VI
 ‡2000-xx



250
 POUND
 WEIGHT CAPACITY
 EVENLY DISTRIBUTED



2010
SERIES
 Wall Saver Arm Chair
 shown in violet | VI
 ‡2010-xx

2000
SERIES
 Wall Saver Arm Chair
 Slant Arms
 shown in violet | VI
 ‡2020-xx



‡ xx in the product number is the placeholder for vinyl/finish choice
 § representative image shows this option ★ options available at additional cost
 1. maximum weight capacity - evenly distributed - refer to user manual 2. accessories sold separately