

MEDICAL SEATING

PHYSICIAN EXAM STOOL

9000 SERIES

Lumex® offers a variety of medical stools to meet the needs of an active facility. Manufactured with the highest quality materials to ensure a long life, Lumex stools offer the most wanted features in an exam stool.

Features:

- Seat diameter comfortably fits to most users
- Pneumatic assist height adjustment for ease of use
- Choice of 3 base styles to complement facility decor
- Choice of height adjustment lever options
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate
- Wide diameter base with hard wheel casters for greater stability
- Backrest, seat height, and footring options
- Available in 24 industry standard colors to complement most decors

Product Info 9501 9702 9903 Fabric Choices



9501 Options:

black nylon composite base +

	Standard Height	with Backrest	with Footring	Lab Height
9501-xx	■			
9501-B-xx	■	■		
9501-FR-xx	■		■	
9501-BFR-xx	■	■	■	
9501-10BFR-xx		■	■	■

★ options available at additional cost



9702 Options:

brushed aluminum base and toe caps +

	Standard Height	with Backrest	with Footring	Lab Height
9702-xx	■			
9702-B-xx	■	■		
9702-FR-xx	■		■	
9702-BFR-xx	■	■	■	
9702-10BFR-xx		■	■	■

★ options available at additional cost



9903 Options:

polished aluminum base +

	Standard Height	with Footring	Lab Height
9903-xx	■		
9903-FR-xx	■	■	
9903-10FR-xx		■	■

★ options available at additional cost

Specifications:

Series	Seat Diameter		Base Diameter with casters		Standard Seat Height				Lab Seat Height				Overall Height w/ Backrest				Footring Diameter		Caster Diameter hard wheel		Shipping Weight ⁸		Weight Capacity ¹ EVENLY DISTRIBUTED		Pneumatic Height Adjustment	Single Lever D-Handle	360° Chrome Handle
	inch	cm	inch	cm	low	high	low	high	low	high	low	high	low	high	low	high	inch	cm	inch	cm	lb	kg ¹	lb	kg ¹			
9501					18.50	46.99	24.00	60.96	24.50	62.23	35.00	88.90	33.00	83.82	39.00	99.06					20	9.07			■	■	
9702	16.00	40.64	24.00	60.96	17.50	44.45	23.00	58.42	23.50	59.69	34.00	86.36	32.50	82.55	37.50	95.25	18.00	45.72	2.00	5.08	22	9.98	350	158.76	■		■
9903					18.00 ²	45.72 ²	23.00 ²	58.42 ²	22.00	55.88	32.00	81.28	18.00 ²	45.72 ²	23.00 ²	58.42 ²					23	10.43			■		■

★ options available at additional cost † -xx in the product number is the placeholder for fabric choice ‡ standard rounding method used
¹ maximum weight capacity - evenly distributed - refer to user manual ² overall height without backrest ⁸ approximate measurement



GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

•1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.

www.grahamfield.com



© 2022 GF Health Products, Inc. All Rights Reserved
 Lumex is a registered trademark of GF Health Products, Inc.

graham field