



Yes, you can.®

6-Point Hydraulic Patient Lift

Model No. 9805P

Safe, Confident
Patient Transfers



Everyday Transfer
Support



Versatile
Low Profile

450 LB WEIGHT
CAPACITY

Versatile Transfers, Multiple Slings

The Invacare 6-Point Hydraulic Patient Lift helps caregivers securely move patients between beds, wheelchairs, recliners, commodes, and other seating surfaces while reducing the physical demands of everyday transfers for caregivers.



Angled Mast Helps
Reduce Patient Sway
During Movement

Adjustable Base Legs
for Improved Access

No Batteries Required

Lightweight
68 Lb Design

Smooth Positioning



The padded 360° 6-point swivel bar rotates to help position patients during transfers while supporting compatibility with multiple sling styles.

Compatible with a Full Portfolio of Invacare Slings

Full body slings • Divided leg slings • Toileting slings

Built to Reduce Caregiver Strain

Thoughtfully designed features help caregivers maneuver and position the lift throughout everyday care routines.

Dual-Side Pump Handle Operation



Rear Shifter Handle for Easy Leg Adjustment

Low-Profile Adjustable Base Legs Open Around Furniture

Low-Friction Wheels Help Reduce Straining



Reliable Daily Support

Built for routine patient transfers in homecare environments with practical features that help caregivers operate the lift with confidence.

- Homecare rehabilitation
- Floor-to-Bed Transfers
- Bed to wheelchair transfers
- Recliner transfers
- Commode transfers





Product Specifications

# Model No. 9805P	5" Low-Friction Swivel Casters
450 lb Weight Capacity	Optional Locking Casters
Manual Hydraulic Lift System	Dual-Side Pump Handle Operation
Padded 360° 6-Point Swivel Bar	Approx. 20" to 64" Lift Range
Adjustable Base Legs	Product Weight: 68 lbs

Tested to meet ISO 10535 standards for strength, durability and safety.

Invacare America
One Invacare Way
Elyria, Ohio 44035-4190
www.invacare.com
(800) 333-6900

©2026 Invacare America. All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare® America or its subsidiaries unless otherwise noted.
Specifications are subject to change without notification. Form No. 26-001, Job. 260104



Yes, you can.®