

CARE CLINER 6530/6540 (w/ NYLON CASTERS)
STANDARD X-LARGE
6531/6541 (w/ STEEL CASTERS)



ISO 9001-2008 Certified

Outstanding Value and Quality in Our Standard Care Cliner

Lower price does not mean lower quality when you buy a Winco Care Cliner. Our standard Care Cliner comes with all the superior features you might find on higher-priced clinical chairs.



KEY FEATURES

- Standard head-rest cover
- Padded armrests
- Dual fold down side-tables
- Removable side panels
- Trendelenburg Positioning



PRODUCT SPECIFICATIONS

CARE CLINER 6530/6540 (w/ NYLON CASTERS) STANDARD X-LARGE 6531/6541 (w/ STEEL CASTERS)

- 1 Standard Trendelenberg release on back with easy-access handle.
- 2 Wide elastic bands on back frame & seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Recessed hand-grips assist in entering and exiting the chair.
- 5 Total-lock and Directional-lock casters provide ease of maneuverability. (Available in Nylon or Steel)
- 6 Removable side panels for easy cleaning and maintenance.



6531
Care Cliner

SPECIFICATIONS

	6530/6531	6540/6541
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	49" (124.5cm)	49" (124.5cm)
Overall Width*	32.25" (81.9cm)*	36.25" (90.1cm)*
Back Height Above Seat	31.5" (80cm)	31.5" (80cm)
Overall Length (Upright)	39" (99cm)	44" (111.8cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	27" (68.6cm)	27" (68.6cm)
Weight	94 lbs. (42.6kg)	111 lbs. (50.4kg)
Weight Capacity	275 lbs. (124.7kg)	450 lbs. (204kg)

*Add 9.25" to chair width for each raised side table

OPTIONS



TB-133 **TB**



LiquiCell **LC**



Heat **HT**
Heat & Massage **HM**



IV Pole **IV**
(Left rear only)



Leg Extensions
RE00



Side Cushion
SC00