



MEDCARE MOBILITY

EMPOWERING INDEPENDENCE. DELIVERING CARE.



www.medcaremobility.com



585-505-6665

HOSPITAL BED MATTRESS COMPARISON CHART

Group A - Basic / Entry-Level Support

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Invacare Deluxe Innerspring Hospital Bed Mattress <i>Firm innerspring</i>	Traditional innerspring hospital mattress	36"	80"	6.5"	350 lbs	Non-powered	Basic	Low	No bariatric option listed.	Users who prefer firm traditional support, basic hospital bed use, and durable everyday mattress support.	View Product
ProBasics Firm Support Innerspring Mattress <i>Firm coil support</i>	Firm support innerspring mattress with Bonnell coil system	36"	80", 84"	6"	300 lbs	Non-powered	Basic	Low	No bariatric option listed.	Homecare, nursing homes, hospitals, and users needing firm restorative support with durable coil construction.	View Product
Protekt Fiber Mattress <i>Firm fiber-core</i>	Fiber-core pressure redistribution mattress	35", 42"	76", 80", 84"	6"	35" models: 350 lbs 42" model: 600 lbs	Non-powered	Basic+	Low	42" wide model available with higher 600 lb capacity.	Firm support needs, lightweight handling, basic pressure redistribution, and facilities needing durable non-powered support.	View Product
Protekt 100 Pressure Redistribution Foam Mattress <i>Soft foam comfort</i>	Soft pressure redistribution foam mattress	36"	76", 80", 84"	6", 9" Plus	350 lbs	Non-powered	Moderate	Basic	No bariatric option listed.	Low-to-moderate pressure risk, elderly homecare users, soft foam comfort, and basic pressure sore prevention.	View Product
Invacare Essentials Mattress <i>Convolved foam</i>	Convolved foam pressure redistribution mattress	36", 42"	80", 84"	Approx. 6"	350 lbs	Non-powered	Moderate	Basic	42" width option available.	Homecare, residential care, elderly users, long-term care, and affordable pressure redistribution support.	View Product

Group B - Better Foam / Prevention Support

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Protekt 300 Pressure Foam Hospital Bed Mattress <i>Multi-tier foam</i>	Pressure redistribution therapeutic foam mattress	36"	76", 80", 84"	6", 9" Plus	350 lbs	Non-powered	High	Moderate	No bariatric option listed.	Moderate-to-high pressure risk, elderly users, rehab patients, heel pressure support, and long-term bed comfort.	View Product
Protekt 500 Gel Foam Medical Mattress <i>Gel foam cooling</i>	VGPT gel-infused pressure redistribution foam mattress	36", 42", 48", 54"	76", 80", 84"	6", 9" Plus	500 lbs	Non-powered	High	Good	Wider therapeutic configurations available up to 54".	Heat-sensitive patients, sweating concerns, long-term bed users, moderate-to-high pressure risk, and advanced foam comfort.	View Product
Invacare Solace Prevention Mattress <i>Three-zone foam</i>	Pressure redistribution therapeutic foam mattress	36"	80", 84"	6", 9" raised side rail option	350 lbs	Non-powered	High	Moderate	No bariatric option on this model.	Pressure sore prevention, heel pressure, body weight distribution, elderly care, and long-term bed support.	View Product

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Joerns PrevaMatt Hospital Bed Mattress Series <i>Multi-model series</i>	Therapeutic pressure redistribution foam mattress series	35", 42"	80", 84"	6", 8" raised perimeter models	Basic: 300 lbs Defend / Console / Alleviate: 500 lbs	Non-powered	Moderate-High	Moderate	42" width options available on select models.	Facilities and homecare buyers needing Basic, Defend, Console, Alleviate, TEMPUR®, heel slope, perimeter, or raised perimeter options.	View Product

Group C - Bariatric Foam Support

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Protekt 600 Bariatric Hospital Bed Mattress <i>Bariatric foam</i>	Bariatric pressure redistribution foam mattress	42", 48", 54", 60"	76", 80", 84", 60" model listed at 80"	6", 9" Plus 60" option listed as 6"	42": 600 lbs 48": 750 lbs 54": 1000 lbs 60": verify by configuration	Non-powered	High	Moderate	Yes - wide bariatric sizes up to 60".	Bariatric pressure redistribution, wider bed users, high body weight support, long-term care durability, and medium-to-high pressure risk.	View Product
Invacare Bariatric Foam Mattress <i>750 lb support</i>	Bariatric high-resilience foam mattress	42"	80"	Standard bariatric therapeutic profile	750 lbs	Non-powered	Moderate	Low-Moderate	Yes - 42" bariatric foam model.	Simple bariatric homecare, heavy-duty everyday support, firm foam comfort, flippable core, and long-term durability.	View Product
Invacare Solace Prevention Bariatric Mattress <i>Built-in side support</i>	Bariatric pressure redistribution foam mattress	42"	80"	Approx. 6"	450 lbs	Non-powered	High	Moderate	Yes - 42" bariatric prevention model.	Bariatric pressure buildup, transfer stability, heel pressure, pressure sore prevention, and stable side support.	View Product
Invacare Softform Premier Bariatric Pressure Ulcer Mattress <i>Premium bariatric foam</i>	Premium bariatric pressure ulcer prevention foam mattress	42", 48"	Standard therapeutic lengths	Premium therapeutic height profile	650 lbs	Non-powered	Very High	Moderate	Yes - 42" and 48" bariatric widths.	Very-high-risk patients, bariatric users, premium foam immersion, pressure ulcer prevention, and transfer edge stability.	View Product

Group D - Air / Alternating Pressure / Low Air Loss

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Protekt Supreme Support Self-Adjusting Air/Foam Mattress <i>Hybrid air/foam</i>	Hybrid self-adjusting air/foam mattress with optional powered pump	35", 42", 48", 54"	80", 84"	7"	35": 500 lbs 42": 600 lbs 48": 650 lbs 54": 700 lbs	Non-powered or powered with optional pump	High-Very High	Good	Yes - wide 42", 48", and 54" options.	Hybrid foam and air support, bariatric users, heel pressure, heat comfort, limited movement, and upgrade flexibility.	View Product

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Protekt Aire 3000 Alternating Low Air Loss Mattress System <i>Accessible air therapy</i>	Alternating pressure and low air loss mattress system	36"	80"	8" 8"/11" with foam raised rails	350 lbs	Powered	High	Good	No bariatric option listed.	Moderate pressure sore risk, moisture concerns, powered therapy, low-pressure awareness, and long-term bed rest.	View Product
Protekt Aire 4000DX Hospital Air Mattress <i>Digital low air loss</i>	Low air loss and alternating pressure mattress system with digital pump	36"	80"	8" 8"/11" with foam raised rails	400 lbs	Powered	Very High	Very Good	No bariatric option listed.	Pressure ulcer treatment and prevention, sweating, moisture buildup, skin maceration, and long-term immobility.	View Product
Invacare microAIR MA500 Alternating Pressure Low Air Loss Mattress <i>Entry powered air</i>	Alternating pressure low air loss mattress system	36"	80"	8"	350 lbs	Powered	High	Good	No bariatric option listed.	Low-to-medium risk patients needing entry-level powered alternating pressure, low air loss, heel slope, and alarms.	View Product
Invacare microAIR MA600 Low Air Loss Mattress <i>3-in-1 therapy</i>	3-in-1 alternating pressure low air loss mattress system	36"	80"	8"	500 lbs	Powered	Very High	Very Good	Higher 500 lb capacity, but no wider bariatric width listed.	Medium-to-high risk patients, heel protection, transfer support, moisture control, and long-term immobility.	View Product
Invacare microAIR MA800 Alternating Pressure Low Air Loss Mattress System <i>True low air loss</i>	Advanced true low air loss, alternating pressure, and pulsation mattress system	36", 42", 48"	80", 84" 42" x 84" listed	10"	36": 600 lbs 42"/48": 1000 lbs	Powered	Very High	Maximum	Yes - 42" and 48" bariatric widths with 1000 lb capacity.	High-risk pressure injury patients, bariatric users, severe moisture management, maceration, sweating, and advanced wound care.	View Product

Group E - Lateral Rotation / Advanced Therapy

Product	Mattress Type	Widths	Lengths	Height	Weight Capacity	Powered?	Pressure Relief	Moisture Control	Bariatric Options	Best For	Link
Protekt Aire 7000 Lateral Rotation Mattress System <i>Rotation + pulsation</i>	Advanced lateral rotation, low air loss, alternating pressure, and pulsation mattress system	36", 42"	80"	8"	660 lbs	Powered	Maximum	Maximum	42" width and 660 lb capacity support larger users.	Pulmonary support, lateral rotation therapy, very-high-risk pressure injury care, severe immobility, and advanced repositioning needs.	View Product
Invacare microAIR MA900 Lateral Rotation Low Air Loss Mattress System <i>Advanced rotation</i>	Lateral rotation true low air loss mattress system	36", 42", 48"	80"	10"	36": 600 lbs 42"/48": 1000 lbs	Powered	Maximum	Maximum	Yes - 42" and 48" bariatric widths with 1000 lb capacity.	Very-high-risk immobile patients, pulmonary complication prevention, bariatric rotation therapy, severe wound prevention, and advanced moisture control.	View Product