



MEDCARE MOBILITY

EMPOWERING INDEPENDENCE. DELIVERING CARE.



www.medcaremobility.com



585-505-6665

PATIENT LIFT SLING COMPARISON CHART

Commode, Toileting & Access Slings

Hygiene access, commode opening, access-style transfer options

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
Invacare Full Body Sling with Commode Opening, Mesh	Full body commode sling	Toileting, shower chair transfers, bed-to-chair, partially or totally dependent users, limited head control, reclined amputee use	Commode opening, quick-drying mesh, washable, designed for easier hygiene transfers	Mesh; capacity not specified	Invacare Reliant 450/600, Jasmine, Berdie Evo XPlus, 9805P lifts
Invacare Mesh Comfort Sling with Commode Opening	Full body / comfort commode sling	Amputees, chronic postural pain, toileting, bathing, users with limited head control	Full head, neck, back, and partial under-thigh support; clear labeling; quick-dry breathable mesh; comfort-focused support	Net / mesh polyester; 500 lb	Invacare Reliant 450/600, Berdie Evo XPlus, Jasmine, 9805P; also noted as compatible with Reliant and I-Lift via 4-point hook-up
McKesson Full Body Commode Sling 4 or 6 Point	Full body commode sling	Toileting, positioning in wheelchair or commode, high-capacity floor-lift use	Commode opening, high back, padded under-thigh section, outer positioning handle, designed to support during transfer procedures	Mesh polyester; 600 lb	Can be used with any floor-style lift; 4 sling points for 4- or 6-cradle point compatibility
Proactive Medical Full Body Mesh Sling with Commode Opening	Full body commode sling	Toileting, showering, bathing, bariatric handling, homecare and facility use	Integrated commode opening, breathable mesh, reinforced stitching, made for hygiene-focused full-body transfers	Polyester mesh; 600 lb	Fits most standard patient lifts
Joerns Hoyer Long Seat Sling with Commode Opening	Long seat commode sling	Paralysis, amputee patients, longer seated use, sitting or lying transfers, toileting access	Long hammock-style support from mid-thigh to upper back, breathable polyester, commode opening, built for extended comfort	Breathable polyester; capacity not specified	Advance, Elara, HPL500, HPL700, Presence, Stature; intended for use with loop spreader bar; 4- or 6-point cradle/spreader bar noted
Joerns Hoyer Comfort Access Padded Sling	Access / padded toileting sling	Toileting and basic transfers, patients who need padded comfort, with or without extra head support depending on option	Adjustable body belt, breathable mesh, extra padding around legs, non-slip back, designed for dressing, undressing and hygiene access	Mesh with padding; capacity not specified	Four clip connections; specifically designed for Joerns Presence, Advance and Stature lifts
Invacare Toileting Sling	Toileting sling	Users with full head and neck control needing padded toileting support in bed or chair situations	Fully padded back and under-thigh comfort, brushed lining, color-coded positioning straps, wet or dry environment use	Polyester fabric; capacity not specified	Invacare Reliant 450/600, Jasmine, 9805P lifts
Joerns Hoyer Professional Access Toileting-Hygiene Loop Style Sling	Access / toileting-hygiene sling	Toileting and hygiene-related passive transfers with loop-style lift support	Access-style design for hygiene positioning; available with or without head support depending on option	Material not specified; 500 lb	Intended for use with a loop-style spreader bar

Full Body, Full Back & General Passive Slings

General passive transfer slings, full support, broader-use models

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
Proactive Medical Full Body Mesh Sling	Full body mesh sling	Home care, bathing, shower support, bariatric handling, wet environments, hospital and rehab transfers	Breathable mesh, wet and dry use, reinforced stitching, machine washable up to 185°F, broad daily transfer use	Polyester mesh; 600 lb	Universally compatible with most floor and ceiling patient lifts
Invacare Mesh Full Body Sling	Full body mesh sling	Bed to wheelchair, toilet, shower chair, floor-to-bed transfers, amputee reclined use, limited head control	Brushed lining to prevent sliding and protect skin, quick-drying mesh, washable, made for dependent or partially dependent users	Mesh; 450 lb	Invacare Reliant 450/600, Jasmine, 9805P lifts
Joerns Hoyer One Piece Sling with Positioning Strap	General-purpose one-piece sling	Basic transfers, hygiene, bathing, amputee support, general homecare and facility use	Positioning strap, soft durable polyester, classic all-in-one design, certified safety build	Polyester; capacity not specified	Designed for Hoyer Classic series and 4-point cradle lifts including HPL450; described as broadly compatible with Hoyer lifts

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
<u>Joerns Hoyer Full Back Sling</u>	Full back sling with head support	Passive lifts, broad general fit, patients needing full back and head support	Anatomically contoured shape, integrated head support, soft leg padding, fits 85–90% of patients	Material not specified; 500 lb	Compatible with 6-point spreader bar lifts like Joerns Hoyer Presence, Stature and Advance
<u>Joerns Hoyer Quickfit Deluxe Sling</u>	General-purpose deluxe sling	General transfers across a very wide patient range	Quickfit-based design for 95% of clients, more supportive lift, available in padded and multiple fabric options, pediatric through XL sizes	Various fabrics; capacity not specified	Used with 6-point cradle on Presence, Stature and Advance lifts

Divided Leg, U-Slings & Universal High Support

More open lower-body positioning, toileting access, high-support variants

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
<u>Invacare Divided Leg Sling</u>	Divided leg sling	Full or partial dependence, toileting access, wheelchair, recliner and shower chair transfers, wet or dry use	Divided leg design, brushed anti-slip lining, full-body support for many users, large toileting opening, padded comfort	Solid composite polyester; 450 lb	Invacare Reliant, Jasmine, Birdie Evo, 9805P lifts; also described as working with Invacare lifts and 6-point spreader bars
<u>Proactive Medical Solid Padded Divided Leg Sling with Head Support</u>	Divided leg sling with head support	Bariatric and elderly transfers, limited neck control, bed-to-chair, toilet and bath transfers, home and facility use	Built-in head support, soft padded fabric, divided leg design, durable for wet or dry use, positioning and total-lift support	Solid padded polyester; 600 lb	Universal design described as compatible with most patient lift systems
<u>Invacare Universal High Sling</u>	Divided leg / universal high sling	Totally dependent, partially dependent, non-weight-bearing, or limited head-control users needing bed, wheelchair, recliner, shower chair, or floor-to-bed transfers	Large toileting opening, full head and neck support in seated and supine positions, adjustable fit, cushioned spacer fabric option for sensitive skin, extra padding and clear attachment points	Premium polyester or spacer fabric; 500 lb	Compatible with Jasmine, Reliant 350, Reliant 450, Reliant 600, and 9805P lifts
<u>Joerns Hoyer Padded U-Sling</u>	Universal padded U-sling with head support	Wheelchair, bed, floor and other general transfers; hygiene and toileting access with a more open lower-body shape	Integrated head support, open U-shape for hygiene access, comfortable secure design, washable to 185°F	Material not specified; capacity not specified	General-use universal design language is used, but exact lift family compatibility is not specified here

Sit-to-Stand & Transfer Slings

Standing support, active participation, shorter transfers

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
<u>Invacare Premier Series Standing Sling for Seniors</u>	Standing / stand-assist sling	Cooperative users with minimal weight-bearing ability and head/neck control who need help transitioning from sitting to standing	Clear labeling, reinforced back support, soft padded leg straps, adjustable fit, easy cleaning	Polyester fabric; capacity not specified	Standing sling category; exact lift family is not specifically named in this listing
<u>Invacare Transport Sling</u>	Transfer / transfer stand-assist sling	Patients who can bear some weight, short-distance transfers, toileting, seated-to-standing transitions	Compact support, brushed lining, dignified toileting and pivot transfers, solid fabric construction	Solid polyester; 440 lb	For use with Invacare Stand-Up lifts; compatible with stand-assist devices like Reliant 350
<u>Joerns Hoyer Deluxe Standing Sling</u>	Sit-to-stand sling	Patients with some weight-bearing ability who need standing transfer support	Non-slip back pad, adjustable waist belt, easy-fit two-point style, fast caregiver application	Material not specified; 440 lb	Designed for Hoyer Journey and Elevate stand-assist lifts
<u>Invacare Premier Series Stand-Assist Sling</u>	Standing / stand-assist sling	Partial mobility, head and neck control, toileting, weight-bearing exercises, quick standing transfers	Reinforced back, labeled attachment points, structured support, quick-fit design	Premium polyester; 450 lb	Designed for use with Invacare Stand-Up lifts
<u>Proactive Medical Sit to Stand Sling</u>	Sit-to-stand sling	Patients who can bear partial weight, rehab, homecare and long-term care standing support	Solid back support, ergonomic attachment, multiple sizes, designed for standing assistance	Nylon / polyester; 600 lb	Described as attaching to most patient lifts and stand-assist devices

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
Invacare Sit to Stand Sling	Sit-to-stand sling	Users who can sit on edge of bed, bend hips, knees and ankles, and bear at least 60% of body weight	Fully padded fabric, soft-brushed lining, stretch-resistant design, color-coded positioning straps	Solid polyester; 440 lb	Invacare Reliant standing patient lift sling; Large and XL size models noted

Amputee, Long Seat & Bariatric Slings

Specialty support for amputee positioning, extended seated use, higher-capacity needs

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
Invacare Comfort Sling for Amputees	Amputee comfort sling	Amputees, horizontal-to-seated transfers, bed to wheelchair, recliner, shower chair, or floor-to-bed support	Designed for amputees, four easy attachment points, adjustable fit, easy maintenance, breathable net fabric	Net fabric; capacity not specified	Four easy attachment points are specifically mentioned; exact lift family is not named in this listing
Joerns Hoyer Longseat Sling	Longseat sling	Amputees and patients who must remain in sling longer during bed making or bathing	Extended-use design, certified safety construction, special-purpose support option	Material not specified; capacity not specified	6-point connection; used with 6-point cradle on Presence, Stature and Advance lifts
Joerns Hoyer Quickfit Bariatric Sling for Calibre and HPL700 Lifts	Bariatric general-purpose sling	Higher-capacity and broader patient needs across a wide range of transfers	Padding around buttocks and leg areas, straight divided leg sections, positioning handles, colored loops, dignity modesty loop, optional padded leg sections	Material not specified; capacity not specified	Fits 6-point cradles and spreader bars; named for Calibre and HPL700 lifts
Joerns Comfort Amputee Sling	Amputee sling	Lower-limb amputees, extended lifting, bathing, bed changes, prolonged support scenarios	Tailored leg positioning, pressure distribution, clip-style stability, professional-grade durability	Breathable polyester; 500 lb	Clip-style attachment; compatible with Joerns 4-point spreader bar systems such as Presence and Stature lifts

Bath, Repositioning, Limb & Disposable Specialty Slings

Specialized products for turning, bathing, limb support, and infection-control workflows

Product	Type	Best Used For	Main Features / What It Is Made For	Material / Capacity	Compatibility / Lift Fit
Joerns Hoyer Limb Sling with Multiple Loop Straps	Limb sling	Targeted limb positioning during patient handling tasks	Multiple loop attachment straps for positioning adjustments; intended for targeted support rather than full-body transfer	Material not specified; 400 lb	Intended for use with a loop-style spreader bar
Joerns Hoyer Repositioning Sling	Repositioning sling	Repositioning, turning, transferring, under-sheet use, extended patient handling support	Breathable pliable mesh, standard loop or figure-of-8 loop options, designed for re-positioning and turning support	Mesh; 1000 lb	Loop style sling; exact lift family is not specifically named
Joerns Hoyer Twin Turner Sling	Turning / repositioning sling	Turning, log rolling, holding patients in position for hygiene, skin care and treatment access	Easy placement and removal, enhanced patient access, designed to simplify caregiver routines	Material and capacity not specified	Exact lift family compatibility is not specifically named
Joerns Hoyer Nylon Mesh Bath Sling	4-point bath sling	Seated bath transfers, wet environments, quick-drying hygiene-related moves	Quick-drying nylon mesh, color-coded sizes, centerpoint marker, with and without head support options noted	Nylon mesh; 600 lb	4-point connection; compatible with HPL700/HPL700WSC, HPL402, P-C-HLA-2, C-HLA-2, C-HLA-2T, C-HLA, C-CB-L2, HML400
Joerns Hoyer Disposable Loop Style Sling	Disposable full back loop-style sling	Single-patient passive lift use, infection control, hospitals and long-term care settings	Contoured easy-fit design, integrated boned head support, leg padding, bulk case format for clean workflow	Disposable / single-patient-use design; 500 lb	Loop-style passive lift sling; exact lift family is not specifically named