



MEDCARE MOBILITY

EMPOWERING INDEPENDENCE. DELIVERING CARE.



www.medcaremobility.com



585-505-6665

PATIENT LIFTING DEVICES COMPARISON CHART

Product	Type	Capacity	Best Use	Transfer Types	Home / Facility Fit	Key Features	Price Range
<u>Protekt 500 Electric Hoyer Lift</u> <i>Heavy-duty electric</i>	Electric full-body Hoyer lift	500 lb	Heavy-duty homecare, bariatric support, floor recovery, low-bed access, and demanding daily transfers.	Bed, wheelchair, chair, low bed, floor recovery, and full-body sling transfers.	Strong fit for homes, hospitals, nursing homes, assisted living, rehab, and facilities needing higher capacity.	6-point spreader bar, floor lift capability, low-bed access, foot pedal base opening, electric controls, emergency lowering, optional digital scale.	Around \$1,255
<u>Protekt Take-A-Long Folding Electric Patient Lift</u> <i>Portable electric</i>	Folding electric patient lift	400 lb	Homecare, travel, compact storage, mobile caregivers, and families needing powered lifting in smaller spaces.	Bed, chair, wheelchair, room-to-room, travel, and full-body sling transfers.	Best for home use, travel use, compact homes, and multi-location caregiving; lighter than many heavy-duty models.	Folding design, battery-powered lifting, tool-free setup, low-battery and overload alerts, battery gauge, emergency stop, optional scale.	Around \$1,079
<u>Protekt Transformer Homecare Patient Lift</u> <i>Manual or upgradeable</i>	Hydraulic patient lift with optional electric upgrade	400 lb	Budget-conscious homecare buyers who want reliable sling-based transfers with the option to upgrade later.	Bed, wheelchair, chair, toilet/commode area, and full-body sling transfers.	Excellent homecare fit; also useful for hospice, rehab recovery, and lower-volume care settings.	Manual hydraulic operation, optional electric upgrade kit, 6-point spreader bar, compact frame, tool-free setup, 3-year frame warranty.	Around \$639 hydraulic / \$999 electric upgrade
<u>Protekt All-In-One Portable Electric Patient Lift</u> <i>Car + floor transfers</i>	Portable electric transfer lift	350 lb	Versatile home or travel transfers where car, floor, chair, wheelchair, and bed movement are major concerns.	Car, floor, bed, wheelchair, chair, and full-body sling transfers.	Strong for homecare, travel, appointments, and caregivers needing more transfer flexibility outside the bedroom.	Foldable frame, built-in battery, 4-point clip-style spreader bar, foot pedal base adjustment, battery alert, manual emergency lowering.	Around \$1,200
<u>Joerns Hoyer HPL402 Patient Lift</u> <i>Premium electric</i>	Electric patient Hoyer lift	400 lb	Smooth electric transfers for homecare, hospitals, rehab, assisted living, and professional caregiving.	Bed, wheelchair, commode, seating, and general full-body sling transfers.	Good fit for home and clinical settings needing trusted electric operation and ergonomic handling.	Adjustable U-base, reverse-mounted actuator for knee clearance, 6-point non-sway cradle, emergency stop, manual emergency lowering, optional scale.	Around \$1,598
<u>Joerns Advance Portable Hoyer Lift</u> <i>Compact folding Hoyer</i>	Portable folding Hoyer lift	341 lb	Portability, compact storage, tight rooms, travel, and transfers where maneuverability matters more than maximum capacity.	Floor, bed, chair, wheelchair, commode-area, and portable transfer situations.	Best for homes, clinics, compact settings, and mobile use; not the top choice for bariatric needs.	No-tools folding design, aluminum frame, swan-neck tapered legs, ergonomic push handle, wide leg opening, EN ISO 10535 tested.	Around \$2,518
<u>Invacare Birdie Evo XPLUS Hoyer Lift</u> <i>Premium facility lift</i>	Premium electric floor lift	463 lb	Hospitals, rehab centers, long-term care, assisted living, and premium homecare needing stable, smooth transfers.	Bed, chair, wheelchair, floor, high surfaces, and full-body sling transfers.	Excellent facility fit; also suitable for homecare buyers wanting a higher-end Invacare lift.	Slow'R boom dampener, ergonomic sling hooks, curved boom for patient eye contact, folding frame, service monitoring, 2 batteries included.	Around \$3,150
<u>Invacare Reliant 600 Heavy-Duty Power Lift</u> <i>600 lb bariatric</i>	Heavy-duty electric bariatric lift	600 lb	Bariatric care, heavier patients, hospitals, long-term care, rehab centers, and high-support transfer routines.	Bed, chair, wheelchair, floor, high surfaces, commode-area, and bariatric sling transfers.	Best for bariatric homecare and professional care facilities requiring heavy-duty lifting capacity.	Extra-wide padded 360° swivel bar, six-point sling hookup, total pinch-point protection, dual 24V batteries, anti-entrapment, manual lowering.	Around \$3,185
<u>Invacare Reliant 450 Battery-Powered Transfer Lift</u> <i>Reliable powered lift</i>	Battery-powered electric patient lift	450 lb	Reliable powered transfers for homecare, facilities, floor lifts, commode access, and repeated caregiver use.	Chair, commode, floor, bed, wheelchair, and high-surface transfers.	Good fit for homes, hospitals, rehab, assisted living, and long-term care settings.	Total pinch-point protection, padded 360° swivel bar, two 24V batteries, manual emergency lowering, anti-entrapment, 150–300 lifts per charge.	Around \$2,180
<u>Joerns Hoyer HPL700 Bariatric</u>	Premium electric bariatric Hoyer	700 lb	High-capacity bariatric care, facilities, hospitals,	Bariatric bed, chair, wheelchair, floor, and	Best for professional bariatric programs and	Powered base, interchangeable cradle options, Hoyer Smart Monitor, optional	Around \$4,626

Product	Type	Capacity	Best Use	Transfer Types	Home / Facility Fit	Key Features	Price Range
Hoyer Lift <i>700 lb bariatric</i>	lift		rehab, long-term care, and larger patient populations.	full-body sling transfers.	facilities needing maximum capacity and advanced monitoring.	integrated scale, ergonomic push handle, EN ISO 10535 tested.	standard / \$7,240 with scale
Joerns Hoyer Presence Medical Hoyer Lift <i>Advanced commercial</i>	Commercial electric medical Hoyer lift	500 lb	High-use care settings needing broad transfer range, powered base, monitoring, and advanced facility performance.	Floor, bed, high surfaces, chair, wheelchair, stretcher option, and full-body sling transfers.	Excellent for hospitals, long-term care, rehab, nursing homes, and advanced homecare setups.	Smart Monitor, powered leg adjustment, interchangeable cradle options, optional scale, optional low leg kit, stretcher option, aluminum frame.	Around \$5,293
Protekt Onyx Hydraulic Patient Lift <i>Budget 450 lb manual</i>	Manual hydraulic patient lift	450 lb	Affordable homecare, elderly transfers, basic facility use, and caregivers wanting a strong manual lift.	Bed, wheelchair, toilet/commode area, recliner, chair, and full-body sling transfers.	Strong homecare and light facility fit; excellent for buyers wanting value and higher manual capacity.	Hydraulic pump, 4- and 6-point sling compatibility, tool-free assembly, hand-operated base, rear locking casters, welded steel frame.	Around \$549
Joerns Hoyer HML400 Manual Hydraulic Lift <i>Classic manual Hoyer</i>	Manual hydraulic Hoyer lift	400 lb	Trusted homecare hydraulic lifting, basic clinical transfers, and caregivers wanting a proven Joerns Hoyer model.	Bed, wheelchair, commode, chair, furniture access, and sling-based transfers.	Best for homecare and light clinical use where manual operation is acceptable.	Rotating pump handle, adjustable U-base, non-sway cradle, disassembles for storage, powder-coated frame, compact homecare design.	Around \$695
McKesson Hydraulic Patient Floor Lift <i>Durable clinical value</i>	Manual hydraulic floor lift	450 lb	Durable homecare or clinical lifting with strong frame construction and affordable manual operation.	Bed, wheelchair, chair, patient floor lift use, and sling-based transfers.	Good for homecare, clinics, facilities, and buyers wanting McKesson durability at a value price.	Silver vein steel frame, 6-point cradle, rotating pump handle, adjustable base, offset mast and boom, 5-inch casters, latex-free, HCPCS E0630.	Around \$649
Invacare Painted Hydraulic Manual Hoyer Lift <i>Affordable Invacare hydraulic</i>	Manual hydraulic Hoyer lift	450 lb	Affordable manual medical lift for homecare, facilities, floor pickup needs, and Invacare sling compatibility.	Floor, bed, chair, wheelchair, commode-area, and full-body sling transfers.	Good homecare and facility fit for buyers wanting an Invacare manual lift at a lower price point.	Six-point padded swivel bar, 360° patient rotation, low-friction 5-inch casters, side-to-side pump handle, wide lift range, disassembles for transport.	Around \$639