



MEDCARE MOBILITY

EMPOWERING INDEPENDENCE. DELIVERING CARE.



www.medcaremobility.com



585-505-6665

WALKERS FOR SENIORS COMPARISON CHART

Standard 4-Wheel Walkers with Seats

Product	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Frame	Best Use
Protekt Bariatric Rollator (500 lb)	500 lbs	24 lbs	18" × 13.5"	22"	34" – 38"	8"	Reinforced steel	Heavy-duty support for bariatric seniors.
Karman R-4608 Rollator	250 lbs	21 lbs (15 lbs frame)	18" × 13.5"	21"	33" – 37"	8"	Aluminum	Lightweight rollator with basket for daily outings.
ProBasics Deluxe Aluminum Rollator	300 lbs	18.4 lbs	14" × 12"	23"	32" – 37.5"	8"	Aluminum	Sturdy lightweight rollator with padded seat & pouch.
ProBasics Aluminum Rollator	300 lbs	15.4 lbs	14" × 12"	21"	30" – 36"	6"	Aluminum	Affordable everyday rollator for seniors.
Carex Classics Steel Rollator	300 lbs	18.5 lbs	12" padded seat	21"	31" – 36"	8"	Steel	Durable budget rollator with bag & locking brakes.
Carex Rollator Steel	350 lbs	17 lbs	12" padded seat	20"	30" – 34.25"	6"	Steel	Steel frame rollator with strong weight capacity.
Carex Four Wheel Walker	300 lbs	17.5 lbs	Padded seat	21"	33.75" – 36.5"	8"	Aluminum	Basic 4-wheel walker with oversized wheels.
Protekt Aluminum Rollator	300 lbs	13 lbs	13.5" × 12"	20"	31" – 36.5"	6"	Aluminum	Ultra-light rollator for indoor/outdoor use.
Protekt Steel Rollator	350 lbs	17 lbs	13.5" × 12"	20"	31" – 36.5"	6"	Steel	Durable steel rollator with flip-up seat.
Lumex Walkabout Imperial Bariatric	450 lbs	~20 lbs	20" × 13.8"	21.8"	34" – 39"	8"	Reinforced aluminum	Wide bariatric rollator with reflectors & pouch.
Lumex Set n' Go Adjustable	300 lbs	~17–18 lbs	14" × 12"	18" – 22" adjustable	29" – 38"	6"	Aluminum	Adjustable height rollator for versatility.
Karman R-4700W Bariatric Steel Rollator	400 lbs	24 lbs	18" seat pad	22"	33" – 37"	8"	Steel	Extra-wide seat & handles for bariatric users.
Karman R-4800 Bariatric Aluminum Rollator	400 lbs	20 lbs	18" seat pad	22"	35" – 39"	8"	Aluminum	Bariatric rollator with compact fold & basket.
Lumex Walkabout Hemi Rollator	300 lbs	14 lbs	14" × 13.5"	20"	32" – 36.5"	6"	Aluminum	Petite-friendly rollator with lower seat height.
Lumex Set n' Go Wide Adjustable	350 lbs	~17–18 lbs	15.75" × 12"	19.5" – 22.5" adjustable	32" – 38"	8"	Aluminum	Wide adjustable seat rollator for extra comfort.

Heavy-Duty & Bariatric Walkers with Seats – Maximum Strength & Support

Product	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Frame Material	Best Use
Protekt Bariatric Rollator	500 lbs	24 lbs	18" × 13.5"	22"	34" – 38"	8"	Steel	Bariatric users needing heavy-duty support indoors & outdoors.
Comodità Prima Rollator Heavy Duty	400 lbs	19.5 lbs	15" × 14.25"	22"	35" – 39"	8"	Reinforced Aluminum	Larger adults needing higher capacity with ergonomic design.
Comodità Prima Special Wide Seat	400 lbs	19.5 lbs	16" × 14.25"	22.5"	35" – 39"	8"	Reinforced Aluminum	Users needing extra-wide orthopedic seat & heavy-duty safety.
Triumph Escape HD Rollator	450 lbs	20.28 lbs	20.67" W	21"	31" – 37.9"	8"	Aluminum	Heavy-duty yet lightweight walker with chair for daily mobility.
Lumex Walkabout Imperial Bariatric Rollator	450 lbs	~20 lbs	20" × 13.8"	21.8"	34" – 39"	8"	Reinforced Aluminum	Bariatric users needing wide seat & reinforced frame.
Karman R-4700W Extra Wide Steel Rollator	400 lbs	24 lbs	18" W seat pad	22"	33" – 37"	8"	Steel	Extra-wide heavy-duty rollator for bariatric daily use.
Karman R-4800 Extra Wide Aluminum Rollator	400 lbs	20 lbs	22" W seat	22"	35" – 39"	8"	Aluminum	Bariatric users wanting lightweight frame & wide comfort seat.

Lightweight & Compact Rollators – Easy to Lift, Fold & Transport

Product	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Frame Material	Best Use
Triumph Essentials Compact Rollator	300 lbs	16.95 lbs	17" × 8"	20"	32" – 36.75"	8" front / 7" rear	Aluminum	Lightweight travel-ready rollator with shopping bag. Ideal for tight indoor spaces and transport.
Triumph Escape HD Rollator	450 lbs	20.28 lbs	20.67" wide seat	21"	31" – 37.9"	8" anti-shock wheels	Aluminum	Heavy-duty yet lightweight rollator with ergonomic seat, storage, and superior capacity.
Comodità Prima Heavy-Duty Rollator	400 lbs	19.5 lbs	15" × 14.25"	22"	35" – 39"	8" wheels	Reinforced Aluminum	Premium heavy-duty rollator with reinforced frame, "S" back design, and

Product	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Frame Material	Best Use
								extra safety for larger users.
Comodità Uno Classic Rollator	300 lbs	19.5 lbs	16" × 13.5"	22"	34" – 39.5"	10" front / 8" rear	Aluminum	Premium orthopedic rollator with wide seat, cup holder, cane holder, and smooth maneuverability indoors & outdoors.
Comodità Spazio Special – Heavy Duty Walker with Seat	440 lbs	20.5 lbs	20" × 14.25"	22"	35" – 39"	8" wheels	Reinforced Steel	Extra-wide orthopedic seat with coccyx relief. Best for bariatric users who need maximum comfort and stability.
COMODITÀ Tipo Classic – Travel Walker with Big Wheels	300 lbs	18.5 lbs	17.5" × 11"	22.5"	31.5" – 37.5"	10" front / 8" rear	Aluminum	Lightweight double-folding design ideal for cruises, air travel, casinos, and everyday outings.

All-Terrain Rollators – Built for Outdoor Adventures

Product	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height Range	Wheel Size	Frame Material	Best Use
COMODITÀ Tipo All Terrain Walker	300 lbs	22 lbs	17.5" × 11"	22.5"	34" – 39"	12" front / 10" rear pneumatic	Reinforced Aluminum	Active users who want an all-terrain outdoor walker with premium accessories.
UNO All Terrain Rollator	300 lbs	28 lbs	16" × 13.5"	22"	34" – 39.5"	12" front / 10" rear pneumatic	Reinforced Aluminum	Rugged outdoor use with orthopedic seat for long walks & adventures.
Comodità Tipo StandUp Advance Upright Walker	300 lbs	29 lbs	17.5" × 11"	21"	41" – 47"	12" front / 10" rear pneumatic	Reinforced Aluminum	Upright posture walking, relieving back/shoulder strain on outdoor terrain.
Lumex Gaitster Upright Rollator	300 lbs	23 lbs	18" × 11"	24"	40.2" – 47"	9.25" front / 8" rear	Aluminum	For seniors needing upright support indoors/outdoors with ergonomic forearm rests.
COMODITÀ Tipo Petite All Terrain	300 lbs	24 lbs	17" × 11"	19"	30" – 34"	10" front / 8" rear pneumatic	Reinforced Aluminum	Petite users up to 5'4" needing outdoor mobility with all-terrain wheels.

Petite & Small User Rollators

Product	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Best Use
Triumph Escape Petite	300 lbs	15.4 lbs	14.5" × 9.5"	19"	28" – 36"	8"	Compact rollator designed for users under 5'4".
COMODITÀ Tipo Petite All Terrain	300 lbs	24 lbs	17" × 11"	19"	30" – 34"	10" front / 8" rear	Petite all-terrain rollator for outdoors & uneven ground.
COMODITÀ Tipo Petite Rollator	300 lbs	18.5 lbs	17.5" × 11"	19"	29.5" – 32.5"	10" front / 8" rear	Lightweight petite walker for tight indoor/outdoor spaces.
COMODITÀ Piccola Petite Rollator	400 lbs	19.5 lbs	15" × 14.25"	18"	29" – 33"	8"	Heavy-duty petite rollator for users up to 5'4", compact yet strong.
Escape Mini Small Rollator	275 lbs	15.6 lbs	18.5" × 9.5"	21"	29.5" – 37.5"	6"	Small frame, lightweight design perfect for petite seniors.

Standard, Shopping & 3-Wheel Walkers Compared

Product	Type	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Best Use	Key Features
Comodità Avanti Walker	Standard Walker	300 lbs	18 lbs	–	–	32" – 38"	5" front wheels	Maximum stability for rehab & post-surgery recovery.	Lightweight reinforced frame, ergonomic grips, stability-focused design.
Triumph Essentials Walker	Standard Walker	300 lbs	17 lbs	–	–	32" – 38"	5" front wheels	Safe indoor mobility, home use, and recovery patients.	Compact fold, lightweight steel frame, ergonomic grips.
Triumph Shopping Rollator	Shopping Rollator	300 lbs	19 lbs	14" × 9.5"	21"	37" – 39.75"	8" front / 7" rear	Perfect for malls, grocery trips & assisted living communities.	Large shopping bag, rear basket, padded seat & backrest, compact fold.
Essentials Shopping Rollator	Shopping Rollator	300 lbs	18 lbs	14" × 9.5"	21"	37" – 39.75"	8" front / 7" rear	Great for seniors needing mobility plus shopping storage.	Large front shopping bag (33 lbs), rear basket, padded seat.
Stander Rollator Walker	Specialty Walker	300 lbs	19 lbs	–	–	Adjustable	8"	Convenient storage-focused design for indoors/outdoors.	Padded seat, ergonomic brakes, built-in storage, compact fold.
Lumex GF 3-	3-Wheel Walker	250 lbs	14 lbs	–	–	33" – 37"	8"	Agile and portable	Zippered pouch,

Product	Type	Weight Capacity	Product Weight	Seat (W × D)	Seat Height	Handle Height	Wheel Size	Best Use	Key Features
Wheel Cruiser								for tight indoor spaces.	foldable frame, ergonomic brakes.
Karman R-3600 3-Wheel Rollator	3-Wheel Walker	250 lbs	15 lbs	–	–	35" – 41"	8"	Great for seniors needing maneuverability & quick folding.	Basket, pouch, tray, quick 3-sec fold, push-down brakes.