



MEDCARE MOBILITY
EMPOWERING INDEPENDENCE. DELIVERING CARE.



www.medcaremobility.com



585-505-6665

WALKERS FOR SENIORS COMPARISON CHART

Product	Walker Type	Description	Frame	Product Weight	Weight Capacity	Seat	Wheel Size	Indoor Use	Outdoor Use	Best For
Group A - Standard Walkers										
Protekt Two-Button Folding Steel Walker	Standard walker	Standard folding medical walker; available with or without wheels	1" steel frame; junior version aluminum	Adult approx. 8-9 lbs; junior with wheels approx. 5 lbs	350 lbs	No seat	5" wheels on wheeled models	Excellent	Limited	Post-surgery, rehab, high stability, walker for surgery, hospital walker, controlled indoor walking.
Protekt Bariatric Two-Button Folding Walker	Bariatric standard walker	Heavy-duty bariatric folding walker with or without wheels	Durable steel frame	12 lbs without wheels; 13 lbs with wheels	500 lbs	No seat	5" wheels on wheeled model	Excellent	Limited	Heavy-duty stability, larger users, bariatric recovery, safer controlled support.
McKesson Bariatric Folding Walker	500 lb folding walker	Bariatric folding walker with 5" wheels and rear glides	Aluminum A-frame with steel legs and side braces	Not listed	500 lbs	No seat	5" front wheels	Excellent	Limited	Bariatric standard walker, larger patients, controlled indoor medical walking support.
McKesson Folding Walker Adjustable Height Aluminum Frame	Adult folding walker	Adult dual-release folding walker with wheels	Aluminum frame	Not listed	350 lbs	No seat	5" front wheels with rear glides	Excellent	Limited	Standard walkers for seniors, post-hospital recovery, indoor mobility, walker patient support.
McKesson Junior Folding Walker	Junior / petite walker	Junior dual-release folding walker with wheels	Aluminum frame	Not listed	350 lbs	No seat	5" front wheels with rear glides	Excellent	Limited	Shorter users, petite elderly, smaller walker size, controlled recovery support.
Group B - Lightweight Rollators / Everyday Walkers With Seat										
Protekt Aluminum Rollator	Lightweight rollator	Lightweight walker for elderly with seat	Lightweight aluminum	13 lbs	300 lbs	13.5" W x 12" D; 20" seat height	6" soft-grip casters	Excellent	Good	Lightweight walkers for seniors, fatigue, shopping, daily use, indoor/outdoor errands.
Protekt Steel Rollator	Steel rollator	Lightweight rollator walker with seat	Durable steel	17 lbs	350 lbs	13.5" W x 12" D; 20" seat height	6" soft-grip casters	Excellent	Good	Affordable senior walker with seat, stronger daily rollator, elderly users needing rest breaks.
Protekt Steel Rollator KD	KD steel rollator	Knock-down mobility walker with seat	Powder-coated steel	17.5 lbs	350 lbs	13.5" W x 12" D; 20" seat height	6" soft-grip casters	Excellent	Good	Reliable mobility walker, compact delivery, tool-free setup, senior daily support.
McKesson 4 Wheel Rollator	Adjustable rollator	Adjustable-height 4-wheel rollator	Aluminum frame	14 lbs	300 lbs	13" seat width; 18"-22" adjustable seat height	6" casters	Excellent	Good	Lightweight walker with wheels, adjustable seat height, daily errands, walkers for handicapped.
Triumph Essentials Adult Walker With Seat	Compact adult walker	Compact adult walker with padded seat	Lightweight aluminum	Not listed	300 lbs	Padded seat with supportive backrest	5" front wheels	Excellent	Light outdoor	Compact adult walker with seat, indoor use, travel storage, simple daily support.
Triumph Essential Compact Lightweight Rollator	Compact rollator	Compact lightweight rollator walker with seat	Lightweight rollator frame	16.95 lbs	300 lbs	17" W x 8" D; 20" seat height	8" front; 7" rear	Excellent	Good	Compact walkers for seniors, travel, car trunk, folding storage, daily errands.
Essential Seated Rollator	Foldable seated rollator	Lightweight foldable rollator walker with seat	Lightweight rollator frame	15.25 lbs	300 lbs	15.5" W x 13.5" D; 22" or 23" seat height	6" wheels	Excellent	Good	Lightweight foldable rollator, daily mobility, travel, senior walker with seat.
Karman R-4608 Rollator	4 wheel walker with seat	Lightweight 4 wheel walker with seat	Powder-coated aluminum	15 lbs frame only; 21 lbs with basket	250 lbs	18" W x 13.5" D; 21" seat height	8" non-marking tires	Excellent	Very good	Best 4 wheel walker with seat, smooth indoor/outdoor wheels, shopping, everyday mobility.
Group C - Shopping, Daily Independence & Stylish Walker Aids										
Essentials Shopping Rollator	Shopping rollator	Adult wheeled walker with seat and shopping storage	Aluminum	18 lbs	300 lbs	14" W x 9.5" D; 21" seat height	8" front; 7" rear	Excellent	Good	Grocery shopping, errands, storage, active seniors, adult wheeled walker with seat.
COMODITA Avanti Walker Aid	Modern walker aid	Modern walker aid for elderly independence	Lightweight steel, foldable	15 lbs	300 lbs	No traditional seat; tray and bag included	8" wheels	Excellent	Good	Stylish walkers for seniors, one-hand braking, stroke weakness, daily home use.
Lumex Set n' Go Wide Rollator	Wide adjustable rollator	Wide adjustable-height 4 wheeler mobility walker	Lightweight aluminum	Approx. 17-18 lbs	350 lbs	15.75" W x 12" D; 19.5"-22.5" adjustable seat height	8" casters	Excellent	Very good	Wider comfort, adjustable fit, daily independence, larger adults, indoor/outdoor use.
Stander Let's Move Rollator	Stylish lightweight rollator	Stylish lightweight rollator	Reinforced aluminum	14.6 lbs	350 lbs	17" W x 7" D faux-leather seat	8" wheels	Excellent	Good	Stylish walkers for adults, lightweight travel, social outings, compact storage.
Group D - Petite / Short User Walkers										
Comodita Piccola Petite Rollator	Petite rollator	Petite rollator walker with seat	Reinforced frame with S-shaped back support	19.5 lbs	400 lbs	15" W x 14.25" D; 18" seat height	8" wheels	Excellent	Good	Petite walkers for elderly, short seniors, small walkers for seniors with seat.
Lumex Walkabout Hemi Rollator	Hemi / petite rollator	Hemi walker with wheels for shorter users	Lightweight aluminum	14 lbs	300 lbs	14" W x 13.5" D; 20" seat height	6" wheels	Excellent	Good	Walker for short seniors, petite adults, lower seat height, lightweight daily use.

Product	Walker Type	Description	Frame	Product Weight	Weight Capacity	Seat	Wheel Size	Indoor Use	Outdoor Use	Best For
COMODITA Tipo Petite Rollator	Petite travel rollator	Petite travel rollator with seat	Reinforced aluminum	18.5 lbs	300 lbs	17.5" W x 11" D; 19" seat height	10" front; 8" rear	Excellent	Very good	Small rollator, petite walker with seat, stylish petite walker, travel and errands.
COMODITA Tipo Petite All Terrain Rollator	Petite all-terrain	Small all-terrain rollator for petite elderly users	Reinforced aluminum	24 lbs	300 lbs	17" W x 11" D; 19" seat height	10" pneumatic front; 8" rear	Good	Excellent	Petite outdoor walker, small rollator for petite elderly, uneven sidewalks and active use.
Group E - Three-Wheel / Compact Walkers										
Karman R-3600 Three-Wheel Rollator	3-wheel rollator	3 wheeled mobility walker	Durable steel	15 lbs	250 lbs	No full sitting seat; tray, basket, and pouch included	Three 8" solid casters	Excellent	Good	Small walker, narrow spaces, apartments, quick errands, compact maneuverability.
Lumex GF 3-Wheel Cruiser	3-wheel cruiser	Lightweight 3 wheeler mobility walker	Lightweight aluminum	14 lbs	250 lbs	No seat; zippered carry pouch included	8" wheels	Excellent	Good	Compact walker, maneuverability, narrow indoor areas, active seniors.
Group F - Heavy-Duty / Bariatric Walkers										
Protekt Bariatric Rollator	500 lb rollator	Heavy duty bariatric rollator with seat	Reinforced steel	24 lbs	500 lbs	18" W x 13.5" D; 22" seat height	8" wheels	Good	Very good	Best heavy duty walker, bariatric users, 500 lb capacity, wide seat support.
Karman R-4800 Extra Wide Aluminum Rollator	Extra wide aluminum	Extra wide heavy duty aluminum rollator	Powder-coated aluminum	20 lbs with basket	400 lbs	22" seat width; 18" seat pad; 22" seat height	8" non-marking casters	Good	Good	Heavy duty rollator, wide seat, bariatric support, lighter aluminum frame.
Karman R-4700W Extra Wide Steel Rollator	Extra wide steel	Extra wide steel bariatric rollator	Super-strength powder-coated steel	24 lbs with basket	400 lbs	23" overall seat; 18" seat pad; 22" seat height	8" non-marking casters	Good	Good	Heavy duty bariatric rollator, wider frame, stronger steel support.
Lumex Walkabout Imperial Bariatric Rollator	Wide seat bariatric	Heavy duty walker with wide seat	Reinforced lightweight aluminum	Approx. 20 lbs	450 lbs	20" W x 13.8" D; 21.8" seat height	8" indoor/outdoor casters	Good	Good	Heavy duty walkers with seats, wide seat comfort, larger adults, bariatric daily use.
Triumph Escape HD	450 lb lightweight HD	Heavy-duty lightweight rollator walker with chair	Lightweight aluminum T-frame	20.28 lbs	450 lbs	20.67" wide; 21" seat height	8" anti-shock swivel wheels	Good	Very good	Heavy-duty yet portable, travel, larger adults, easy brakes, strong rollator with chair.
Comodita Prima Special	Wide seat rollator	Wide seat walker with 4 wheels and seat	Reinforced frame with S-shaped wheel support	19.5 lbs	400 lbs	16" W x 14.25" D; 22.5" seat height	8" rubber wheels	Good	Good	Wide seat rollator, comfort, safety, indoor/outdoor use, larger adults.
Comodita Spazio Heavy Duty Walker With Seat	440 lb heavy duty	Heavy duty walker with seat	Reinforced heavy-duty steel frame	20.5 lbs	440 lbs	18" W x 14.25" D; 22" seat height	8" rubber tires	Check width	Good	Walkers for heavy people, extra room, strong seating support, 30" doorway recommended.
Comodita Spazio Special	Extra-wide orthopedic	Premium heavy duty walker with extra-wide orthopedic seat	Reinforced heavy-duty frame	20.5 lbs	440 lbs	20" W x 14.25" D; 22" seat height	8" rubber tires	Needs 30" doorway	Good	Extra-wide comfort, bariatric users, orthopedic seat support, heavy-duty stability.
Group G - Outdoor / All-Terrain / Big Wheel Walkers										
COMODITA Tipo All Terrain Walker	All-terrain walker	Premium all-terrain rollator walker	All-terrain rollator frame	22 lbs	300 lbs	17.5" W x 11" D; 22.5" seat height	12" pneumatic front; 10" pneumatic rear	Good	Excellent	Outdoor walkers for elderly, grass, gravel, sidewalks, trails, parks, active seniors.
COMODITA Tipo Classic Travel Walker	Big wheel travel walker	Travel walker with big wheels	Patented compact travel rollator frame	18.5 lbs	300 lbs	17.5" W x 11" D; 22.5" seat height	10" front; 8" rear	Good	Very good	Lightweight travel walker, big wheels, car trunks, sidewalks, travel, active users.
UNO All Terrain Rollator	Premium outdoor rollator	Premium outdoor all-terrain rollator	Heavy-duty reinforced frame with S-shaped rear support	28 lbs	300 lbs	16" W x 13.5" D; 22" seat height	12" pneumatic front; 10" pneumatic rear	Good but larger	Excellent	Outdoor mobility, gravel, trails, boardwalks, active seniors, rugged all-terrain use.
Stander Let's Go Out Rollator	Outdoor stylish rollator	Outdoor-friendly stylish rollator	Lightweight rollator frame	17.45 lbs	300 lbs	17" x 7"; 23.6" seat height	10" front; 8" rear with suspension	Good	Very good	Stylish outdoor rollator, travel, errands, smoother outdoor surfaces, social use.
Group H - Upright Walkers & Hybrid Rollator / Transport Chair										
Comodita Tipo StandUp Advance	Upright all-terrain	Upright walker with seat	Premium upright all-terrain frame	29 lbs	300 lbs	17.5" W x 11" D; 21" seat height	12" pneumatic front; 10" pneumatic rear	Good but larger	Excellent	Upright posture, forearm support, wrist/back relief, outdoor travel, large wheels.
Lumex Gaitster Upright Rollator	Upright rollator	Upright rollator walker with seat	Upright rollator frame	23 lbs	300 lbs	18" W x 11" D; 24" seat height	9.25" front; 8" rear	Good	Very good	Upright posture support, reduced wrist strain, taller walking position, forearm support.
Karman R-4602-T Two-in-One Rollator and Transport Chair	Rollator transport chair	Rollator and transport chair hybrid	Rollator transport frame	21 lbs listed; 28 lbs shipping/spec listing	250 lbs	13" seat width; 22" seat height	8" flat-free casters	Good	Good	Airports, travel, caregiver pushing, fatigue, rollator with footrest, long outings.