

BARIATRIC HOME BED



FULL ELECTRIC ULTRA WIDE 48" and 54" BED FRAME

Reinforced steel tubular deck for extra support and durability



1,000 LBS.

**WEIGHT
CAPACITY**

Color Coded Actuator Control Box

Easily locate where to
attach the actuators and
hand pendant



3 Linear Actuators for Head, Foot and Hi/Lo functions

Sturdy 2 piece reinforced tubular deck

Illuminated 8 function hand pendant

4 x 5" castors - 2 locking and 2 swivel

Durable and washable bed ends

Bed Dimensions: 88" L x 48" or 54" W x 18" - 27 H

Side rails - half rail



ILLUMINATED 8 FUNCTION HAND PENDANT

- Raise and lower the head deck
- Raise and lower the foot deck
- Raise and lower the height of the bed
- Auto contour function

LOW POSITION: 18"

HIGH POSITION: 27"

HEAD ANGLE: 0 - 60 Degrees

FOOT ANGLE: 0 - 35 Degrees

WEIGHT CAPACITY: 1,000 Lbs.

MODEL #: HCFE48-HR

MODEL #: HCFE54-HR