

MATTRESS OVERLAY



BUBBLE PAD

ALTERNATING PRESSURE PUMP AND PAD SYSTEM

**DESIGNED TO BE USED IN THE
PREVENTION, TREATMENT AND
MANAGEMENT OF PRESSURE ULCERS**



Adjustable
Pressure



PRESSURE RANGE:
50 - 105mmHg



Whisper Quiet
Pump

CYCLE TIME: 10 mins.

AIR OUTPUT: 4 - 6 liters

PAD SIZE: 79" x 35" x 3"

**130 Individual
Alternating
Bubble Cells**

MODEL # **BP300**

