

WHEELCHAIR BARIATRIC 26"/28"/30" WIDTH



HEAVY DUTY, ULTRA WIDE, BARIATRIC AVAILABLE IN 26", 28" & 30" WIDTH

WEIGHT CAPACITY: 700 LBS.



- **DUAL CROSSBAR FOR EXTRA SUPPORT**
- **DETACHABLE ARM & FOOT REST WITH ELEVATED LEG REST**
- **8" DURABLE FRONT SWIVEL CASTER WITH REINFORCED ALUMINUM FORK**



- Detachable desk length armrest
- Padded armrest
- Elevated leg rest, adjustable height with aluminum footrest
- Detachable foot rest
- Push to lock chrome steel brakes
- Silver vein steel frame with steel side panels
- 24" Composite Rear Mag Wheels
- 8" Heavy Duty Front Swivel Casters
- Aluminum heavy duty reinforced front fork
- Dual crossbar supports heavy weight
- Nylon padded Flame Retardant upholstery

Model #	Seat Dimensions	Seat Height	Back Height	Total Height	Weight Capacity
WC2620	26" W X 20" D	19"	18"	90"	700 LBS.
WC2820	28" W X 20" D	19"	18"	90"	700 LBS.
WC3020	30" W X 20" D	19"	18"	90"	700 LBS