

GEL AND FOAM WHEELCHAIR CUSHION



The **Medacure Premium Foam and Gel Wheelchair Cushions** are designed to redistribute pressure and weight, aiding in the prevention, treatment, and management of pressure ulcers. These cushions provide enhanced support and comfort for individuals who remain seated for long periods. Featuring a waterproof, breathable 2- way stretch cover with an anti-slip bottom, they ensure durability and stability. Available in a range of standard and bariatric sizes to suit diverse needs.

- **FOAM:** High density polyurethane foam
- **GEL:** Offers excellent pressure redistribution and efficiently conducts heat, delivering a cooling sensation during prolonged sitting.



2" Height




3" Height

GEL AND FOAM WHEELCHAIR CUSHION



- Helps in the prevention of pressure sores
- Aerated foam promotes coolness and dryness for the user
- Ensures proper positioning and enhances patient comfort
- Two-way stretch nylon cover offers breathability, waterproof protection, and comfort
- Includes four convenient tie straps to securely attach the cushion to the wheelchair
- Features a replaceable cover for hassle-free cleaning
- Warranty: 1 Year
- Weight Capacity:
Standard: 250 lbs. / Bariatric: 350 lbs.



Replaceable Cover For Hassle-Free Cleaning



Dual-compartment gel bladder prevents gel migration

Item #	Description/Size	Case Pack
GEL CUSHION 2"		
GC1616-2	Gel Cushion 16" x 16" x 2"	4
GC1816-2	Gel Cushion 18" x 16" x 2"	4
GC1818-2	Gel Cushion 18" x 18" x 2"	4
GC2016-2	Gel Cushion 20" x 16" x 2"	4
GC2018-2	Gel Cushion 20" x 18" x 2"	4
GEL CUSHION 3"		
GC1616-3	Gel Cushion 16" x 16" x 3"	4
GC1816-3	Gel Cushion 18" x 16" x 3"	4
GC1818-3	Gel Cushion 18" x 18" x 3"	4
GC2016-3	Gel Cushion 20" x 16" x 3"	4
GC2018-3	Gel Cushion 20" x 18" x 3"	4
GC2216-3	Gel Cushion 22" x 16" x 3"	4
GC2218-3	Gel Cushion 22" x 18" x 3"	4
GC2418-3	Gel Cushion 24" x 18" x 3"	4
GC2618-3	Gel Cushion 26" x 18" x 3"	4
GC2818-3	Gel Cushion 28" x 18" x 3"	4
GC3018-3	Gel Cushion 30" x 18" x 3"	4

Item #	Description/Size	Case Pack
FOAM CUSHION 2"		
FC1616-2	Foam Cushion 16" x 16" x 2"	4
FC1816-2	Foam Cushion 18" x 16" x 2"	4
FC1818-2	Foam Cushion 18" x 18" x 2"	4
FC2016-2	Foam Cushion 20" x 16" x 2"	4
FC2018-2	Foam Cushion 20" x 18" x 2"	4
FOAM CUSHION 3"		
FC1816-3	Foam Cushion 18" x 16" x 3"	4
FC2016-3	Foam Cushion 20" x 16" x 3"	4
FC2018-3	Foam Cushion 20" x 18" x 3"	4
FC2216-3	Foam Cushion 22" x 16" x 3"	4
FC2218-3	Foam Cushion 22" x 18" x 3"	4
FC2418-3	Foam Cushion 24" x 18" x 3"	4
FC2618-3	Foam Cushion 26" x 18" x 3"	4
FC2818-3	Foam Cushion 28" x 18" x 3"	4