

FULL ELECTRIC HOMECARE BED

MODEL #: HCFE36

MODEL #: HCFE36 - LOW


HCPCS CODE: E0297

The Medacure Full Electric Homecare bed is lightweight and easy to install. Available in a standard and low height bed. The 3 Function hand pendant raises the head, foot and height of the bed. The mattress deck has reinforced steel slats for safety and comfort with a built in mattress keeper. Bed ends have a stylish walnut finish. Includes emergency hand crank in case of power failure. Full or Half rails included.



Low Position	15" (Standard Bed) 9" (Low Bed)
High Position	23"
Head Angle	0 - 70 Degrees
Foot Angle	0 - 45 Degrees
Weight Capacity	450 Lbs

- Sturdy reinforced slat deck provides patient comfort and safety
- Built-in emergency hand crank lowers and raises the bed
- 4 x 3" castors – 2 locking and 2 swivel
- Durable and washable bed ends
- Bed Dimensions: 86" L x 36" W
- Mattress retainer keeps mattress from sliding
- Side rails – Full or Half Rail optional



MULTI FUNCTION HAND PENDANT

1. Raise and lower the head deck
2. Raise and lower the foot deck
3. Raise and lower the height of the bed



HAND
CRANK
FUNCTION