

# 3 FUNCTION FULL & SEMI ELECTRIC HOME CARE BED



HALF RAIL



FULL RAIL



Hand Crank for HI/LO function  
on Semi Electric Bed



Mattress retainer holds  
mattress firmly in place

Sturdy 2 piece reinforced tubular deck

Built-in emergency hand crank lowers  
and raises the bed

4 x 3" castors - 2 locking and 2 swivel

Durable and washable bed ends

Bed Dimensions: 88" L x 36" W x 15" - 23 H

Side rails - Full or half rail optional

**LOW POSITION:** 15"

**HIGH POSITION:** 23"

**HEAD ANGLE:** 0 - 70 Degrees

**FOOT ANGLE:** 0 - 45 Degrees

**WEIGHT CAPACITY:** 450 Lbs.



## MULTI FUNCTION PENDANT

- Raise and lower the head deck
- Raise and lower the foot deck
- Raise and lower the height of the bed - Full Electric

MODEL #: HCFE36 - FULL ELECTRIC

MODEL #: HCSE36 - SEMI ELECTRIC