


LINCOLN EXPANDABLE BARIATRIC BED



Bariatric



LINCOLN EXPANDABLE
BARIATRIC BED

TRUE BARIATRIC FRAME

EXPANDS IN THE WIDTH FROM 36" - 42" - 48"
AND IN THE LENGTH FROM 80" - 88"



**LOCKS SECURELY IN
PLACE WITH A HEAVY
DUTY SPRING BOLT**

**BUILT-IN
EXPANSION**

**NO TOOLS
REQUIRED**



**10 FUNCTION HAND
PENDANT WITH
UNIVERSAL SYMBOLS**

Hi/Low, Head,
Foot, Trendelenburg,
Reverse Trend



Expands & Retracts
without Tools!



Heavy duty rectangle
tube bed deck with
reinforced steel frame



10 Position hand pendant
with Trendelenburg Function



Tool-less Design -
Width & Length expand
& retract without tools




8 - 3" Durable Casters:
4 Swivel and
4 Locking Casters




Expands in length from 80" - 88"

- Expandable width from 36" - 42" - 48"
- Expandable length from 80" - 88"
- Expandable headboard from 36" - 42" - 48"
- Pair of Quarter rails or assist bars - included / Head and Foot boards - included
- IV pole mount and Trapeze bracket mount on 2 head corners
- High: 30" / Low: 14" / Also available in an Ultra Low Bed - Low: 9" / High: 25"
- Under bed clearance of 5.7"
- **WEIGHT CAPACITY: 750 Lbs.**

MODEL # **LX-BARI**
MODEL # **LX-BARI-9 (LOW BED)**



Mahogany



Amber Cherry



Light Oak

