

LINCOLN EXPANDABLE BARIATRIC BED



Bariatric

With BUILT-IN DIGITAL SCALE SYSTEM



ULTIMATE PATIENT SAFETY

- Weighs Patient in Bed
- Large Easy to Read Color Digital Display Panel

ILLUMINATED SCALE SYSTEM with Nurse Control & Patient Override Function

- Nurse Control and Pendant control all 10 positions
- Lockout feature on nurse control overrides patient hand pendant
- Large illuminated color screen displays patient weight
- Stores weight history and calculates difference in weight
- LB./KG selectable
- CPR and emergency stop function
- Comes pre-calibrated
- Powered by AC or backup 12V rechargeable battery included

MODEL # **LX-BARI-S**

DIGITAL SCALE / NURSE CONTROL FEATURES & FUNCTIONS



Digital Scale System – weighs patient in the bed



2 white LED's under the nurse control illuminate the foot end at night



Bed alarm function – sounds when patient is off the bed



Battery power indicator on display

LINCOLN EXPANDABLE BARIATRIC BED

WITH SPLIT FRAME DESIGN

Ideal Rental Solution - folds compact for easy transport and storage



MODEL # **LX-BARI-SF**

Model #	Width	Length	Low Range	High Range	Weight Cap
LX-BARI	36" - 42" - 48"	80" - 88"	14"	30"	750
LX-BARI-9 (Ultra Low)	36" - 42" - 48"	80" - 88"	9"	25"	750
LX-BARI-S (Built-In Scale)	36" - 42" - 48"	80" - 88"	14"	30"	750
LX-BARI-SF	36" - 42" - 48"	80" - 88"	14"	30"	750