

LONG TERM CARE BED ULTRA WIDE



ULTRA LOW BED



AVAILABLE IN STANDARD AND
ULTRA LOW



8 Function
Hand Pendant
With Auto
Contour
Function

STANDARD AND ULTRA LOW BED

**ULTRA WIDE
DESIGNED FOR MAXIMUM
PATIENT COMFORT**

42" FIXED WIDTH MATTRESS DECK

Expandable length from 80" - 88"



Quick Release
Head & Foot
boards included



2 Adjustable
Side Rails included




Durable
4" Locking Casters




Weight Capacity:
SLB: 550 Lbs. /
ULB: 450 Lbs.

- One piece reinforced steel frame
- 8 Function patient hand pendant included
- IV pole socket on all 4 bed corners
- Durable rectangle tube bed deck
- Linak motors - dependable & ultra quiet
- **UNDERBED CLEARANCE: 5.7"**



Mahogany



Amber Cherry



Light Oak



Model #	Description	Width	Length	Low Range	High Range	Weight Cap
SLB42	American Spirit	42" Fixed Width	80" - 88"	14"	30"	550
ULB42	Ultra Low	42" Fixed Width	80" - 88"	9"	30"	450