

FLOOR STANDING BARIATRIC TRAPEZE



DESIGNED TO ASSIST PATIENT'S WHILE
CHANGING POSITION IN BED AND AS AN AID FOR
TRANSFERRING IN AND OUT OF BED



Rubber Coated Grab Bar
allows for a Secure Grip



Adjustable Chain Height
Locks Securely in Place



Adjust Chain Position
from Edge of Headboard
to end of Boom



Raised Rear Casters for
Easy Transport



Rubber Coated Safeguards
keeps Base Secure

LIGHT WEIGHT AND EASY TO ASSEMBLE

- **BOOM LENGTH:** 37"
- **INSIDE WIDTH BETWEEN LEGS:** 37.5"
- **OUTSIDE WIDTH BETWEEN LEGS:** 41.5"
- **BASE LEGS LENGTH:** 49"
- **TOTAL HEIGHT:** 74"
- **WEIGHT CAPACITY:** 500 Lbs.

MODEL # **TR500**