

LONG TERM CARE BED SUPER LOW TO FLOOR



TRAVELS ALL THE WAY DOWN TO FLOOR LEVEL

Ultimate Fall Prevention Solution

Available in 36" and 42" Ultra Wide Bed Deck



80" Wide / 20" High
Underbed Clearance



**10 FUNCTION
HAND PENDANT**
Hi/Lo, Head, Foot,
Trendelenburg,
Reverse Trend

TRAVEL RANGE: 3.9" - 25"

PATIENT HANDSET:
10 function with trendelenburg

DURABLE EASY ROLL CASTERS:
3" with brake

MATTRESS RETAINERS
On each corner - included

WEIGHT CAPACITY: 450 Lbs.

MATRESSS DECK: 80" x 36"

PATIENT ASSIST BAR:
Set of 2 - included



ACTUATOR COVER PROTECTOR

Prevents damage to Motors
& Risk of Operation failure
ensures a Safe Environment
for patient and staff



SPLIT FRAME BED DESIGN

Excellent rental solution
for hassle free transport

TOTAL FOLDED DIMENSIONS:
20" W X 36" or 42" L X 48.5" H



3" Locking Casters
All 4 Casters Lock

MODEL # **ULB3.9**
MODEL # **ULB3.9-42**