

LONG TERM CARE BED ULTRA LOW & HIGH

★★★★★ EXPANDABLE WIDTH FROM 36" - 42" - 48"



Valu Bari Bed

3 POSITION EXPANDABLE WIDTH
36" - 48"

EXTENSIVE TRAVEL RANGE:
7" - 30"


EXPANDABLE LENGTH:
80" - 88"

8 Function Hand Pendant

HI/LO, HEAD, FOOT,
AUTO CONTOUR



Expands & Retracts
without tools!



LOW POSITION OF ONLY 7"!
IDEAL SOLUTION FOR FALL PATIENTS



**LOCKS SECURELY IN
PLACE WITH A HEAVY
DUTY SPRING BOLT**

**BUILT-IN
EXPANSION**

**NO TOOLS
REQUIRED**

- 8 Function hand pendant with universal symbols
- Auto contour function
- Expandable head board expands from 36" to 48" to prevent gaps
- Quick Release boards
- Head and Foot Boards included
- Bumper bar protects walls from scratches and scuff marks
- Durable reinforced steel tube frame
- IV sockets on bed corners
- Corner mattress retainers included
- **8 - 3" CASTERS:** 4 locking and 4 swivel
- **ULTRA HEAVY WEIGHT CAPACITY:** 550 Lbs.



Expandable width from 80" - 88"



Quick Release **Expandable**
Head & Foot Boards



Model #	Description	Expandable Width	Length	Low Range	High Range	Weight Cap
ULB48-X-7/30	Ultra Low	36" - 42" - 48"	80" - 88"	7"	30"	550