

LONG TERM CARE BED ULTRA LOW & HIGH



WITH CENTRAL LOCKING SYSTEM 36" FIXED WIDTH

ULTRA LOW BED

VERSATILE PERFORMANCE & TRAVEL RANGE



8 Function Hand Pendant

HI/LO, HEAD, FOOT, AUTO CONTOUR



EXTENSIVE TRAVEL RANGE:
7" - 30"

EXPANDABLE LENGTH:
80" - 84"



LOW POSITION OF ONLY 7"! IDEAL SOLUTION FOR FALL PATIENTS




CENTRAL LOCKING SYSTEM

Color coded, dual pedal floor lock system keeps bed securely in place when engaged

- 8 Function hand pendant with universal symbols
- Auto contour function
- Quick Release Head & Foot boards
- Head and Foot Boards included
- Bumper bar protects walls from scratches and scuff marks
- Durable reinforced steel tube frame
- IV sockets on bed corners
- Corner mattress retainers included
- **8 - 3" CASTERS:** 4 locking and 4 swivel with directional lock
- One touch, dual pedal, color coded floor lock system
- **HEAVY WEIGHT CAPACITY:** 500 Lbs.



Directional steering lock



Quick Release Head & Foot Boards



Model #	Description	Width	Length	Low Range	High Range	Weight Cap
ULB7/30-CLS	Ultra Low with Floor Lock	36"	80" - 84"	7"	30"	500