

WHEELCHAIR FULL RECLINING HIGHBACK



- **PADDED ELEVATED LEG REST**
- **HIGH BACK with NECK SUPPORT CUSHION**
- **ANTI-TIPPER WHEELS**
- **AVAILABLE IN 18", 20", 22" & 24" WIDTHS**

1. Full Reclining high back with Neck Support Cushion
2. Detachable Padded Armrest
3. Detachable, Swing-Away and Elevated Leg rest
4. Deluxe padded vinyl upholstery with heavy duty inner, flame retardant and mildew resistant.
5. 8" Front Casters with solid rubber tires, aluminum front forks adjustable in three positions.
6. 24" Composite Mag wheels
7. Push-to-lock wheel locks
8. Anti-Tipper wheels
9. Steel side panels



MODEL # 1818-R

WHEELCHAIR FULL RECLINING HIGHBACK



Detachable High Back
Converts to Standard
Wheelchair



Dual Handbrakes



Reclining Easily Adjusts
and Locks in Place



Full Position Reclining



Removeable Headrest

Model #	Overall Length	Overall Width	Overall Height	Folded Size	Seat Width	Seat Depth	Back Height	Seat Height To Floor	Weight	Weight Cap
W1816-R	49"	26"	50"	12"	18"	16"	30"	20"	58 LBS.	300 LBS.
W2016-R	49"	27.8"	50"	12"	20"	16"	30"	20"	59 LBS.	300 LBS.
W2218-R	51"	30"	51"	13.5"	22"	18"	32"	20"	65 LBS.	450 LBS.
W2418-R	51"	32"	51"	13.5"	24"	18"	32"	20"	66 LBS.	450 LBS.