

WHEELCHAIR STANDARD WIDTH



LIGHTWEIGHT AND MAINTENANCE FREE
AVAILABLE IN 18" & 20" WIDTH



Dual Axle
Adjustable Seat Height



Large Storage Pocket
on Back of Seat



Removable & Adjustable
Elevated Leg Rest



Detachable Arm & Foot Rest
Available with a Swing Away & Elevated Foot Rest

WHEELCHAIR EXTRA WIDE BARIATRIC



22" & 24" WIDE

DUAL CROSSBAR

500 LB. WEIGHT CAPACITY



8" Durable Front Caster
Swivel & Adjustable



Dual Crossbar for Extra Support

- Deluxe padded vinyl upholstery with heavy duty inner, flame retardant and mildew resistant
- Desk length and full length armrest available
- Padded detachable armrest
- Elevated leg rest, adjustable height
- Swing away and detachable foot rest
- Push to lock brakes
- Powdered coated steel frame with steel side panels
- 24" Composite Rear Mag Wheels
- 8" Front Casters with solid rubber tires, aluminum front forks adjustable to 3 positions
- Dual Axle - adjustable seat height - allows for repositioning for hemi seat-to-floor height
- Easy-lift rear bar aids in lifting over bumps and curb

Item #	Seat Dimensions	Weight Cap
W1816	18" W x 16" D	300 LBS.
W2016	20" W x 16" D Dual Crossbar	350 LBS.
W2218	Bariatric - 22" W x 18" D Dual Crossbar	500 LBS.
WB2418	Bariatric - 24" W X 18" D Dual Crossbar	500 LBS.

WHEELCHAIR FULL RECLINING HIGHBACK

STANDARD & BARIATRIC ★★★★★



- **PADDED ELEVATED LEG REST**
- **HIGH BACK with NECK SUPPORT CUSHION**
- **ANTI-TIPPER WHEELS**
- **AVAILABLE IN 18", 20", 22" & 24" WIDTHS**

1. Full Reclining high back with Neck Support Cushion
2. Detachable Padded Armrest
3. Detachable, Swing-Away and Elevated Leg rest
4. Deluxe padded vinyl upholstery with heavy duty inner, flame retardant and mildew resistant.
5. 8" Front Casters with solid rubber tires, aluminum front forks adjustable in three positions.
6. 24" Composite Mag wheels
7. Push-to-lock wheel locks
8. Anti-Tipper wheels
9. Steel side panels



WHEELCHAIR FULL RECLINING HIGHBACK

STANDARD & BARIATRIC ★★★★★



Detachable High Back
Converts to Standard
Wheelchair

Dual Handbrakes



Reclining Easily Adjusts
and Locks in Place

Full Position Reclining

Removeable Headrest

Model #	Overall Length	Overall Width	Overall Height	Folded Size	Seat Width	Seat Depth	Back Height	Seat Height To Floor	Weight	Weight Cap
W1816-R	49"	26"	50"	12"	18"	16"	30"	20"	58 LBS.	300 LBS.
W2016-R	49"	27.8"	50"	12"	20"	16"	30"	20"	59 LBS.	300 LBS.
W2218-R	51"	30"	51"	13.5"	22"	18"	32"	20"	65 LBS.	450 LBS
W2418-R	51"	32"	51"	13.5"	24"	18"	32"	20	66 LBS.	450 LBS.

WHEELCHAIR BARIATRIC 26"/28"/30" WIDTH



HEAVY DUTY, ULTRA WIDE, BARIATRIC AVAILABLE IN 26", 28" & 30" WIDTH

WEIGHT CAPACITY: 700 LBS.

- **DUAL CROSSBAR FOR EXTRA SUPPORT**
- **DETACHABLE ARM & FOOT REST WITH ELEVATED LEG REST**
- **8" DURABLE FRONT SWIVEL CASTER WITH REINFORCED ALUMINUM FORK**



- Detachable desk length armrest
- Padded armrest
- Elevated leg rest, adjustable height with aluminum footrest
- Detachable foot rest
- Push to lock chrome steel brakes
- Silver vein steel frame with steel side panels
- 24" Composite Rear Mag Wheels
- 8" Heavy Duty Front Swivel Casters
- Aluminum heavy duty reinforced front fork
- Dual crossbar supports heavy weight
- Nylon padded Flame Retardant upholstery

Model #	Seat Dimensions	Seat Height	Back Height	Total Height	Weight Capacity
WC2620	26" W X 20" D	19"	18"	90"	700 LBS.
WC2820	28" W X 20" D	19"	18"	90"	700 LBS.
WC3020	30" W X 20" D	19"	18"	90"	700 LBS